

- + Backed by 800 Published Studies
- + Fully-Featured Physiological Waveforms
- + AccuWin Pro<sup>TM</sup> 4 ABPM Analysis & Reporting
- + Industry-Leading Measurement Performance

## + SphygmoCor Inside

The gold standard in Pulse Wave Analysis & Central Blood Pressure Measurement

Validated to inter-arterial catheter

Over 800 published studies

Sub-Diastolic CBP Waveform Capture

Fully featured central blood pressure waveforms



## + Advanced Analysis with AccuWin Pro™ 4

Automated interpretation of ABPM data

Programmable quality control metrics

Detailed pulse wave analysis (PWA)

HL7 Output Capability

Automated White-coat analysis

#### + SunTech® ABPM Blood Pressure Cuffs

Patented stretch sleeve keeps the cuff in place both day and night

Optimal patient comfort and data validity

4 cuff sizes to fit almost any arm (18cm-55cm)

Removable bladder for machine washing

### + Specifications

**Data Storage** 

Validations BHS (A/A), ESH, ANSI/AAMI SP10:2002

**BP Technique** Oscillometry with step deflation

**BP Range** 25-260 mmHg **HR Range** 40-200 bpm

Sampling Periods3 adjustable (awake, asleep, and special)Time Intervals5, 10, 15, 20, 30, 45, 60, 90 and 120 min.Event KeyProgram to allow additional patient readings

Flash memory stores up to 250 readings

**Power Source** Two AA batteries

Weight 284 grams, including batteries

**Dimensions** 12 x 7 x 3 cm

PC Requirements Serial or USB port, Windows® XP or later

# +The Oscar 2™ with Sphygmocor Inside System Includes

Oscar 2<sup>TM</sup> ABP monitor

AccuWin Pro™ 4 software

2 SunTech® ABPM cuffs (Sizes 2 and 3)

Carrying pouch & patient belt

USB data cable

4 AA alkaline batteries

User's Manual & Quick Start Guide

Part Numbers	Descriptions
99-0133-03	Oscar 2 <sup>™</sup> with SphygmoCor Inside ABPM System
98-0239-01	ABPM Cuff Size 1 (18-26 cm)
98-0239-02	ABPM Cuff Size 2 (25-35 cm)
98-0239-03	ABPM Cuff Size 3 (32-44 cm)
98-0239-04	ABPM Cuff Size 4 (42-55 cm)



SunTech Medical, Inc. 507 Airport Boulevard, Suite 117 Morrisville, NC 27560-8200 USA

Tel: + 1.919.654.2300 1.800.421.8626 Fax: + 1.919.654.2301 SunTech Medical, Ltd.
Oakfield Industrial Estate
Stanton Harcourt Road
Eynsham, Oxfordshire OX29 4TS
England

Tel: + 44 (0) 1865.884.234 Fax: + 44 (0) 1865.884.235

EC REP

www.SunTechMed.com

SunTech Medical, Ltd.

105 HuanGuan South Road, Suite 15 2~3/F

DaHe Community Guanlan, LongHua District, Shenzhen GuangDong PRC 518110

Tel: + 86-755-2958 8810

+ 86-755-29588986 (Sales)

+ 86-755-29588665 (Service)

Fax: + 86-755-2958 8829 www.SunTechMed.com.CN

